



ICU bed

Motorised ICU bed

- Remote Operated adjustment of
 - Hi-Low
 - Back
 - Knee
 - Trendlenburg /Reverse Trendlenburg
- * Frame are made by Crca M.S Tubes and powder coated frame.
- * ABS butterfly Side-Rails With Angle Indicator.
- * 4 Section top, made by Perforated Crca Ms Sheet
- * SS Telescopic IV rod with four hooks and fitted on 4 locations of bed.
- * Urine Bag holder.
- * Movable on 5" medical wheels with 2 wheel diagonally lockable.
- * Overall Approx. Size:- 2100Mm (L) X 1000Mm (W) X 500-800Mm (H)