

floSKAN 3000

Micro-continuous Syringe Pump



Salient Features

- 3.5" Backlit Color LCD Touch screen
- Multiple Infusion modes
- 11 occlusion pressure levels
- Advanced transmission & multi-layered safety to ensure accurate infusion
- Automatic detection of syringes of various sizes
- Compatible with 5, 10, 20, 30, and 50 (60) ml syringes
- Large infusion range, from 0.1 to 2000 ml/hr
- Built-in library of 1000 drugs
- 6 Hour Battery backup, upgradable to 12 hours
- Designed in accordance with IEC and ISO guidelines



RADIOLOGY



MONITORING



CARDIOLOGY



ANAESTHESIA



VENTILATION



DENTAL



SURGERY



TELE SOLUTIONS




VETERINARY



Easy operation and safety

- User friendly interface
- Stable drive mechanism
- Accurate 16-step motor
- Built-in Drug library
- 3 Bolus modes
- Safety handle to prevent accidental infusion
- Recent Therapy administered



Dual CPU control system, better control and data transmission



Unique Multi function port for DC input, nurse call & software upgrade



Independent Audio & Visual alarm system

Specification

Power Supply:	AC- 100-240V; External DC- 12V; built-in lithium battery, at least 6 hours
Compatible syringe:	5ml, 10ml, 20ml, 30m and 50 (60) ml
Infusion rate:	0.1-2000ml/h, 0.01ml/h increment*
Infusion mode:	Rate, Time, Body weight, Sequence & Relay mode
Infusion accuracy:	Mechanism accuracy: $\pm 1\%$, Accuracy with syringe: $\pm 2\%$
Bolus rate:	0.1- 2000ml/h
KVO:	0.1- 5ml/h

* 0.01ml/h increment till 100m/hr, 0.1ml/h increment from 100-1000m/hr and 1ml/h increment from 1000-2000m/hr



CUSTOMER INTERACTION CENTER

Toll Free No.: **1800-425-7002**
Email : cic@skanray.com

Monday to Saturday - 10AM to 6:45PM

Distributor:

Skanray Technologies Limited

Regd Office & Technology Division: #15-17 Hebbal Indl Area, Mysore - 570016, Karnataka, INDIA. Ph: +91 821 2415559, Fax: +91 821 2403344.
Healthcare Division: Plot #360, KIADB Industrial Area, Hebbal, Mysore - 570018, Karnataka, India. Ph: +91 821 2407000, Fax: +91 821 2407001.
Skanray Europe S.R.L: Via Cicogna 36 & 38, 40068 San Lazzaro di Savena, Bologna, ITALY. Phone: +39 051 6255019.
Email: enquiry@skanray.com | www.skanray.com | CIN: U72200KA2007PLC041774

Warning: The equipment is to be used by only qualified Medical Personnel and in accordance with instructions of the manufacturer, as specified in the User Manual. It should be serviced and maintained by the Company's authorized service personnel. Skanray continuously upgrades its technology, performance and safety features. Pictures and specifications featured in this catalog are only indicative. Please refer to the commercial offer document for exact specifications and performance.

